

CERTIFICATE OF ANALYSIS

Feline Panleukopenia Virus (FPLV)

FA Substrate Slide

Catalog No.:	SLD-IFA-FPL
Lot:	P260515
Expiration:	06 July 2034

Quality Control Method: Indirect FA using CPV positive control (catalog no. PC-IFA-CPV-G) and anti-canine IgG affinity purified FITC conjugate (catalog no. CJ-F-CANG-AP-10ML).

Specific Reaction: 3–4+ signal with the positive control and negative with the diluent control. No background. There are 0–12 infected cells per high-power field.

Other Comments: NA

QC Manager:



Rose Lucas

Pattern Of Fluorescence: Single cells with cytoplasmic and/or nuclear fluorescence.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

13 July 2026