

CERTIFICATE OF ANALYSIS

Canine Adenovirus (CAV)

FITC 50X Conjugate

Catalog No.:	CJ-F-CAV-50X
Lot:	P260531
Expiration:	20 July 2034

Quality Control Method: Direct FA on CAV-2 12-well in-house slides.

Specific Reaction: 3-4+ signal with trace background at 1/50 dilution.

Other Comments: NA

QC Manager:



TJ Yockachonis

Pattern Of Fluorescence: Mostly large, single nuclear inclusion and some cytoplasmic fluorescence in a few infected cells.

FA Cross-Reaction Screening: The raw material has also been screened by FA and has been found to react with Canine Adenovirus Type 1 (CAV-1) with 2–4+ signal, no background at 1/16. It does not react with Canine Coronavirus (CCV), Canine Distemper Virus (CDV), Canine Herpesvirus Type 1 (CHV-1), Canine Parainfluenza Virus Type 2 (CPI-2), Canine Parvovirus (CPV), Rabies Recombinant Nucleoprotein (rNP), Reovirus (REO), Toxoplasma gondii, Vesicular Stomatitis Virus Indiana and New Jersey strains (VSV).

Technical Data Sheet Version: 1

FOR *IN VITRO* LABORATORY USE ONLY.

23 July 2026