

CERTIFICATE OF ANALYSIS

## Feline Herpesvirus Type 1 (FHV-1/FVR)

FA Substrate Slide

|              |              |
|--------------|--------------|
| Catalog No.: | SLD-IFA-FVR  |
| Lot:         | P260603      |
| Expiration:  | 02 July 2034 |

**Quality Control Method:** Indirect FA using FHV-1/FVR positive control (catalog no. PC-IFA-FVR), and anti-feline IgG affinity purified FITC conjugate (catalog no. CJ-F-FELG-AP-10ML).

**Specific Reaction:** 4+ signal with the positive control and negative with the diluent control. No background. There are 0–30 infected cells per high-power field.

**Other Comments:** NA

**QC Manager:**



Rose Lucas

**Pattern Of Fluorescence:** Individual rounded cells with cytoplasmic and nuclear fluorescence. Some small syncytia.

Technical Data Sheet Version: 2

FOR *IN VITRO* LABORATORY USE ONLY.

13 July 2026